


1. Part Number: 675 7009A01

2. Photo:



3. Product Information:

Description	Corer
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Corer Overall Length:	230mm
Corer Overall Diameter:	12mm
Corer Internal Diameter:	9mm
Max. Sampling Depth:	205mm
Push Rod Length:	265mm
Nominal Weight - Corer:	150g
Nominal Weight – Push Rod:	190g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: The Corer and the Push Rod are supplied together.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24