

1. Part Number: 675 A105100

2. Photograph:



3. Product Information:

Description	PTFE Funnel 100mm
Material of Construction:	PTFE
Construction:	Single piece, crevice free
Nominal Max. Outer Diameter:	104 mm
Nominal Outer Diameter of Spout:	22 mm
Nominal Length of Spout:	75 mm
Nominal Overall Height:	150 mm
Nominal Weight:	133 g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Aug-25