

1. Part Number: 675 A10550

2. Photograph:



3. Product Information:

Description	PTFE Funnel 50mm
Material of Construction:	PTFE
Construction:	Single piece, crevice free
Nominal Max. Outer Diameter:	52 mm
Nominal Outer Diameter of Spout:	16 mm
Nominal Length of Spout:	50 mm
Nominal Overall Height:	84 mm
Nominal Weight:	23 g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Aug-25