

1. Part Number: 675 8023P250

2. Photograph:



3. Product Information:

Description:

Powder Funnel - 250mm

Material of Construction:

Polypropylene (PP) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Maximum Diameter:

287mm

Overall Height:

200mm

Outside Diameter of Tip of Spout:

55mm

Inside Diameter of Tip of Spout:

50mm

Length of Spout:

30mm

Nominal Weight:

195g

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25