

1. Part Number: 675 8111H200

2. Photograph:



3. Product Information:

Description:	SteriWare Powder Funnel 200mm
Use:	Single use, disposable powder funnel
Material of Construction:	High Density Polyethylene (HDPE) - virgin • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC Regulations 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Diameter:	200mm
Overall Height:	165mm
OD at Bottom of Spout:	43mm
ID at Bottom of Spout:	40mm
Nominal Weight:	86g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A