

1. Part Number: 675 8111H50

2. Photograph:



3. Product Information:

Description:	SteriWare Powder Funnel 50mm
Use:	Single use, disposable powder funnel
Material of Construction:	High Density Polyethylene (HDPE) - virgin • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC Regulations 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Diameter:	68mm
ID at Top of Funnel:	50mm
Overall Height:	45mm
OD at Bottom of Spout:	18mm
ID at Bottom of Spout:	15mm
Nominal Weight:	5g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture