

**1. Part Number:** 675 8370H100

**2. Photograph:**



**3. Product Information:**

Description:

SteriWare Funnel 100mm

Use:

Single use, disposable funnel

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Max. Diameter:

117mm

ID at Top of Funnel:

100mm

Overall Height:

103mm

OD at Bottom of Spout:

11mm

ID at Bottom of Spout:

8mm

OD at Top of Spout:

17mm

Nominal Weight:

16g

Moulding Packing Environment:

Class 100,000 Cleanroom

Primary Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

N/A

Quantity per Box:

50

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested