

1. Part Number: 675 A244250

2. Photograph:



3. Product Information:

Description:	Industrial Liquid Funnel
Material of Construction:	304 Stainless Steel
Top Diameter:	240mm
Outer Diameter of Spout:	20mm
Overall Length:	235mm
Nominal Spout Length:	88mm
Nominal Weight:	330g
Side Handle:	Spot welded to cone section of funnel
Recommended Storage Conditions:	Dry and ambient temperature

Note: Due the manufacturing method there is an internal crevice where the cone and spout are joined. These funnels are primarily designed for industrial applications.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24

1. Part Number: 675 A244310

2. Photograph:



3. Product Information:

Description:	Industrial Liquid Funnel
Material of Construction:	304 Stainless Steel
Top Diameter:	310mm
Outer Diameter of Spout:	20mm
Overall Length:	325mm
Nominal Spout Length:	88mm
Nominal Weight:	530g
Side Handle:	Spot welded to cone section of funnel
Recommended Storage Conditions:	Dry and ambient temperature

Note: Due the manufacturing method there is an internal crevice where the cone and spout are joined. These funnels are primarily designed for industrial applications.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24