

1. Part Number: 675 5100E1000

2. Photo:



3. Product Information:

Description

PTFE Extension Rod 1000mm

Product Contact Parts:

PTFE (Polytetrafluoroethylene)

- Pure, virgin PTFE
- USP Class VI compliant
- FDA conforming

Overall Length:

1000mm

Working Temperature Range:

-200 to +280°C

Nominal Weight:

430g

Recommended Storage Conditions:

Dry and ambient temperature

Notes:

This Extension Rod is suitable for use with any of the following PTFE Dippers:

- 5100A-100
- 5100A-250
- 5100A-500
- 5100A-1000

Multiple Extension Rods can be screwed together to further increase the reach of the sampling cup

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Dec-24