

1. Part Number: 675 7010A02

2. Photo:



3. Product Information:

Description	Maxi Corer 600mm
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.8 microns Ra
Diameter of Outer Tube:	25mm
Overall Length:	600mm
Blade Length:	500mm
Handle Length	100mm
Nominal Weight of Sampler:	600g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24