

1. Part Number: 675 8064H500S

2. Photograph:



3. Product Information:

Description:	Powder Probe 600mm - Sterile Single use
Use:	powder sampler
Material of Construction:	High Density Polyethylene (HDPE) – virgin <ul style="list-style-type: none">• Conforms to FDA 21 CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Ultrasonic welding
Nominal Length:	615mm
Blade Length:	497mm
Diameter of Blade:	22mm
Width of opening in blade:	10mm
Nominal Weight of Sampler:	89g
Moulding Environment:	Class 100,000 Cleanroom
Welding & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested