

1. Part Number: 675 6011A20D

2. Photo:



3. Product Information:

Description	Discharge Rod for Powered Ice Sampler 20mm
Materials:	316L stainless steel
Method of Construction:	Screw together
Surface Finish:	Better than 1.0 microns Ra
Overall Length:	255 mm
Body Diameter:	14 mm
Nominal Weight:	440g
Recommended Storage Conditions:	Dry and ambient temperature

Note: This device is designed to work with the Powered Ice Sampler (Part No. 6011A-20).

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25