

1. Part Number: 675 8008P2000

2. Photograph:



3. Product Information:

Description:	SteriWare Graduated Beaker 2000ml
Use:	Single use, disposable beaker
Material of Construction:	Polypropylene (PP) - virgin • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	2000ml
Brimful Volume:	2400ml
Max. Diameter (inc. spout):	165mm
Top diameter (excl. spout)	146mm
Bottom Diameter:	121mm
Overall Height:	195mm
Graduation Marks:	From 100ml to 2000ml (every 100ml)
Nominal Weight:	161g
Moulding Environment:	Class 100,000 Cleanroom