

1. Part Number: 675 8008P25

2. Photograph:



3. Product Information:

Description:	SteriWare Graduated Beaker 25ml
Use:	Single use, disposable beaker
Material of Construction:	Polypropylene (PP) - virgin • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	25ml
Brimful Volume:	35ml
Max. Width (inc. spout):	43mm
Top diameter (excl. spout)	38mm
Bottom diameter:	33mm
Overall Height:	46mm
Graduation Marks:	From 5ml to 25ml (every 5ml)
Nominal Weight:	5g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	50