

1. Part Number: 675 8008P250

2. Photograph:



3. Product Information:

Description:	SteriWare Graduated Beaker 250ml
Use:	Single use, disposable beaker
Material of Construction:	Polypropylene (PP) - virgin • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	250ml
Brimful Volume:	320ml
Max. Diameter (inc. spout):	80mm
Top diameter (excl. spout)	74mm
Bottom Diameter:	60mm
Overall Height:	108mm
Graduation Marks:	From 50ml to 250ml (every 10ml)
Nominal Weight:	33g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A