

1. Part Number: 675 8008P50

2. Photograph:



3. Product Information:

Description:	SteriWare Graduated Beaker 50ml
Use:	Single use, disposable beaker
Material of Construction:	Polypropylene (PP) - virgin • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	50ml
Brimful Volume:	65ml
Max. Width (inc. spout):	53mm
Top diameter (excl. spout)	48mm
Bottom diameter:	41mm
Overall Height:	54mm
Graduation Marks:	From 10ml to 50ml (every 5ml)
Nominal Weight:	8g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)