

**1. Part Number:** 675 A1791000

**2. Photograph:**




### 3. Product Information:

Description:	Aluminium Bottle 1000ml
Material of Construction:	Bottle    Aluminium Cap       Polypropylene (PP)
Nominal Volume:	1000ml
Brimful Volume:	1200ml
Overall Height:	245mm
Max. Bottle Diameter:	89mm
ID of Bottle Neck:	27mm
Nominal Weight:	143g (bottle & cap)
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Feb-20