

1. Part Number: 675 A179250

2. Photograph:




3. Product Information:

| | |
|---------------------------------|---|
| Description: | Aluminium Bottle 250ml |
| Material of Construction: | Bottle Aluminium Cap Polypropylene (PP) |
| Nominal Volume: | 250ml |
| Brimful Volume: | 290ml |
| Overall Height: | 160mm |
| Max. Bottle Diameter: | 55mm |
| ID of Bottle Neck: | 21mm |
| Nominal Weight: | 48g (bottle & cap) |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Feb-20