



1. Part Number: 675 1030B all lengths

2. Photograph:



3. Product Information:

Description:

PowderThief - 19mm diameter

Part No:	Overall Length	Weight
675 1030B600	680mm	1075g
675 1030B1000	1080mm	1580g
675 1030B1500	1580mm	2080g
675 1030B2000	2080mm	2620g

Material of Construction:

316L stainless steel
(except for the Clutch which is made from PTFE)

External Polish:

<1 microns Ra

Internal Polish:

None

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25