


1. Part Number: 675 1030BT all sizes

2. Product Information:



Part No:	Nominal Volume	Weight
675 1030BT2	2ml	190g
675 1030BT3	3ml	180g
675 1030BT5	5ml	168g
675 1030BT8	8ml	115g
675 1030BT10	10ml	128g

Description

Tips for use with 19mm diameter PowderThief Bodies

Materials of Construction:

316L stainless steel

Finish:

<1 microns Ra

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin

Recommended Storage Conditions:

Dry and ambient temperature

3. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25