



**1. Part Number:** 675 1030A all lengths

**2. Photograph:**



**3. Product Information:**

Description: PowderThief - 12.5mm diameter

Part No:	Overall Length	Weight
675 1030A300	370mm	350g
675 1030A600	665mm	540g
675 1030A1000	1065mm	795g
675 1030A1500	1565mm	1020g
675 1030A2000	2065mm	1485g

Material of Construction: 316L stainless steel (except for the Clutch which is made from PTFE)

External Polish: <1 microns Ra

Internal Polish: None

Recommended Storage Conditions: Dry and ambient temperature

Spare Parts: PTFE Clutch Part No: SS002300  
Clutch Housing Part No: SS001900

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 9-Jan-25