


1. Part Number: 675 1030AT all sizes

Description


Tips for use with 12.5mm diameter PowderThief Bodies

Photo	Part No:	Nominal Volume	Weight
	675 1030AT025	0.25ml	18g
	675 1030AT050	0.50ml	19g
	675 030AT075	0.75ml	20g
	675 1030AT100	1.00ml	24g
	675 1030AT150	1.50ml	26g
	675 1030AT200	2.00ml	31g

Materials of Construction:

Tip
Collar

316L stainless steel
PTFE (Part No. SS002000)



Finish:

<0.5 microns Ra

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are free of animal products

Recommended Storage Conditions:

Dry and ambient temperature

2. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25