

1. Part Number: 675 1210J2500

2. Photo:



3. Product Information:

Description	Aluminium Slot Sampler 2500mm
Material of Construction:	Anodised aluminium
Method of Construction:	Tip riveted to outer tube
Overall Length:	2570mm
Max. Sampling Depth:	1855mm
Number of Slots:	7
Slot Length:	100mm
Max. Slot Width:	27mm
Body Diameter:	50mm
Max. Sample Volume:	2370ml
Nominal Weight:	3.6 Kg
Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 13-Nov-24

1. Part Number: 675 1210K1000

2. Photo:



3. Product Information:

Description: Jumbo Slot Sampler 1000mm

Material of Construction: 316L stainless steel

Outer Diameter of Outer Tube: 40 mm

Inner Diameter of Inner Tube: 32 mm


Overall Length: 1050 mm

Slot Length: 120 mm

Slot Width: 20 mm

Number of Slots: 3

Slot Positions:



External Polish: <1 Micron Ra

Internal Polish: None

BSE/TSE statement: All polishing compounds used in the manufacture of this product are of vegetable origin

Nominal Weight of Sampler: 3.4 Kg

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-24

1. Part Number: 675 1210K1500

2. Photo:



3. Product Information:

Description Jumbo Slot Sampler 1500mm

Material of Construction: 316L stainless steel

Outer Diameter of Outer Tube: 40 mm

Inner Diameter of Inner Tube: 32 mm


Overall Length: 1550 mm

Slot Length: 120 mm

Slot Width: 20 mm

Number of Slots 4

Slot Positions:



External Polish: <1 Micron Ra

Internal Polish: None

BSE/TSE statement: All polishing compounds used in the manufacture of this product are of vegetable origin

Nominal Weight of Sampler: 4.9 Kg

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-24

1. Part Number: 675 1210K2000

2. Photo:



3. Product Information:

Description: Jumbo Slot Sampler 2000mm

Material of Construction: 316L stainless steel

Outer Diameter of Outer Tube: 40 mm

Inner Diameter of Inner Tube: 32 mm


Overall Length: 2060 mm

Slot Length: 120 mm

Slot Width: 20 mm

Number of Slots: 4

Slot Positions:



External Polish: <1 Micron Ra

Internal Polish: None

BSE/TSE statement: All polishing compounds used in the manufacture of this product are of vegetable origin

Nominal Weight of Sampler: 6.5 Kg

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-24