

1. Part Number: 675 1210K1000

2. Photo:



3. Product Information:

Description: Jumbo Slot Sampler 1000mm

Material of Construction: 316L stainless steel

Outer Diameter of Outer Tube: 40 mm

Inner Diameter of Inner Tube: 32 mm


Overall Length: 1050 mm

Slot Length: 120 mm

Slot Width: 20 mm

Number of Slots: 3

Slot Positions:



External Polish: <1 Micron Ra

Internal Polish: None

BSE/TSE statement: All polishing compounds used in the manufacture of this product are of vegetable origin

Nominal Weight of Sampler: 3.4 Kg

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-24