


1. Part Number: 675 1220J1500

2. Photo:



3. Product Information:

Description	Jumbo Pocket Sampler 1500mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Internal Polish:	None
Overall Length:	1570mm
Diameter of Outer Tube:	40mm
Pocket Length:	125mm
Pocket Width:	20mm
Nominal Pocket Volume:	90ml
Number of Pockets:	4
Slot Positions:	



Nominal Weight:	5.6 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-24