

**1. Part Number:** 675 1360A1501

**2. Photo:**



### 3. Product Information:

Description: AccuPod Sample Thief 1500mm, 1 port  
Material of Construction: 316L stainless steel  
Method of Construction: Welds ground and polished  
Surface Finish: Better than 0.8 microns Ra (tube exteriors)  
Overall Length: 1500mm  
Outer Tube Diameter: 25mm  
Number of Sampling Ports: 1  
Position of Sampling Ports:



Nominal Weight: 3.8 Kg  
BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin  
Recommended Storage Conditions: Dry and ambient temperature  
First Use: This sampler should be thoroughly cleaned before being used for the first time  
Note: This sampler is designed for use with AccuPod inserts. AccuPod inserts must be ordered separately. Contact QAQC LAB for details.

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24