


1. Part Number: 675 1360A301

2. Photo:



3. Product Information:

Description:	AccuPod Sample Thief 300mm, 1 port
Material of Construction:	316L stainless steel
Method of Construction:	Welds ground and polished
Surface Finish:	Better than 0.8 microns Ra (tube exteriors)
Overall Length:	300mm
Outer Tube Diameter:	25mm
Number of Sampling Ports:	1
Position of Sampling Ports:	



Nominal Weight:	960g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature
First Use:	This sampler should be thoroughly cleaned before being used for the first time
Note:	This sampler is designed for use with AccuPod inserts. AccuPod inserts must be ordered separately. Contact QAQC LAB for details.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24