


**1. Part Number:** 675 1360A603

**2. Photo:**



### 3. Product Information:

Description:	AccuPod Sample Thief 600mm, 3 ports
Material of Construction:	316L stainless steel
Method of Construction:	Welds ground and polished
Surface Finish:	Better than 0.8 microns Ra (tube exteriors)
Overall Length:	600mm
Outer Tube Diameter:	25mm
Number of Sampling Ports:	3
Position of Sampling Ports:	



Nominal Weight:	1.6 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature
First Use:	This sampler should be thoroughly cleaned before being used for the first time
Note:	This sampler is designed for use with Accupod inserts. AccuPod inserts must be ordered separately. Contact QAQC LAB for details.

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24