


**1. Part Number:** 675 5035J1500

**2. Photograph:**



**3. Product Information:**

Description

Jumbo Liquid Master – 1500mm

Materials of Construction:

The product contact parts are either PTFE or FEP

- Tip PTFE
- Outer Tube FEP
- Stopper PTFE
- Cup PTFE
- Spring Mild steel (non-product contact part)
- Inner Rod PTFE & stainless steel
- Knob PTFE (non-product contact part)

Method of Construction:

The threaded parts of the device screw together. The device can be quickly dis-assembled for cleaning

Overall Length:

1565mm

Diameter of Outer Tube:

35mm

Max. Diameter of Tip:

38mm

Diameter of Hole in Tip:

9mm

Nominal Weight:

1355g

Recommended Storage Conditions:

Dry and ambient temperature

**4. Approval:**

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 2-Dec-24