

1. Part Number: 675 5065B2002

2. Photo:



Angular Beaker fitted to a Telescopic Rod.
Note: Telescopic Rod is sold separately.



Angular Beaker 2000ml
(Blue rotating elbow is included in the price)

3. Product Information:

Description	Angular Beaker 2000ml	
Materials of Construction:	Beaker	Polypropylene - Translucent
	Rotating Elbow	PBT – Blue
Nominal Volume:	2000ml	
Brimful Volume:	2400ml	
Maximum Beaker Diameter:	145mm	
Overall Beaker Width:	176mm	
Nominal Beaker Height:	195mm	
Nominal Weight:	240g	
Recommended Storage Conditions:	Dry and ambient temperature	

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 11-Dec-24