

**1. Part Number:** 675 5065C1001B

**2. Photo:**



Pendulum Beaker fitted to a Telescopic Rod.  
Note: Telescopic Rod is sold separately.



Pendulum Beaker 1000ml  
(Blue rotating elbow is included in the price)

**3. Product Information:**

Description	Pendulum Beaker 1000ml – 316Lss
Materials:	Beaker: 316L stainless steel Rotating Elbow PBT - Blue
Surface Finish:	Better than 0.8 microns Ra
Nominal Volume:	1000ml
Internal Beaker Diameter:	100mm
Overall Beaker Diameter:	110mm
Nominal Beaker Height:	145mm
Nominal Weight:	750g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: Suitable for use with 5065A and 5065S telescopic rods.

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 11-Dec-24