

1. Part Number: 675 5065C1002

2. Photo:



Pendulum Beaker with Telescopic Rod.
Note: Telescopic Rod is sold separately.



Pendulum Beaker
(Beaker + Rotating Elbow)

3. Product Information:

Description

Pendulum Beaker 1000ml

Materials of Construction:

Beaker	Polypropylene - Translucent
T-Bar	Polypropylene - Blue
Rotating Elbow	PBT - Blue

Nominal Volume:

1000ml

Maximum Beaker Diameter:

115mm

Overall Beaker Width:

145mm

Nominal Beaker Height:

148mm

Nominal Weight:

200g

Recommended Storage Conditions:

Dry and ambient temperature

Note: The Rotating Elbow will fit into the shank on the side of the Beaker.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Jun-24