

**1. Part Number:** 675 5065C252

**2. Photo:**



Pendulum Beaker with Telescopic Rod.  
Note: Telescopic Rod is sold separately.



Pendulum Beaker  
(Beaker + Rotating Elbow)

### 3. Product Information:

|                                 |                             |                             |
|---------------------------------|-----------------------------|-----------------------------|
| Description                     | Pendulum Beaker 250ml       |                             |
| Materials of Construction:      | Beaker                      | Polypropylene - Translucent |
|                                 | T-Bar                       | Polypropylene - Blue        |
|                                 | Rotating Elbow              | PBT – Blue                  |
| Nominal Volume:                 | 250ml                       |                             |
| Maximum Beaker Diameter:        | 73mm                        |                             |
| Overall Beaker Width:           | 1010mm                      |                             |
| Nominal Beaker Height:          | 109mm                       |                             |
| Nominal Weight:                 | 132g                        |                             |
| Recommended Storage Conditions: | Dry and ambient temperature |                             |

Note: The Rotating Elbow will fit into the shank on the side of the Beaker.

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Jun-24