

1. Part Number: 675 5065C502

2. Photo:



Pendulum Beaker with Telescopic Rod.
Note: Telescopic Rod is sold separately.




Pendulum Beaker
(Beaker + Rotating Elbow)

3. Product Information:

Description	Pendulum Beaker 500ml	
Materials of Construction:	Beaker	Polypropylene - Translucent
	T-Bar	Polypropylene - Blue
	Rotating Elbow	PBT – Blue
Nominal Volume:	500ml	
Maximum Beaker Diameter:	92mm	
Overall Beaker Width:	120mm	
Nominal Beaker Height:	112mm	
Nominal Weight:	145g	
Recommended Storage Conditions:	Dry and ambient temperature	

Note: The Rotating Elbow will fit into the shank on the side of the Beaker.

4. Approval:

Signed:


Name: Kuldeep Flora
Position: Quality Manager
24-Jun-24