

1. Part Number: 675 5065D502

2. Photo:



3. Product Information:

Description

Bottle Holder & 500ml bottle

Materials of Construction:

Bottle Holder Fame	Polypropylene	Blue
Strap	Nylon 66	Black
Locking Wheel	Acetal	Black
Rotating Elbow	PBT	Blue
Bottle	HDPE	Natural
Bottle Lid	Polypropylene	Natural

Nominal Bottle Volume:

500ml

Maximum Width of Bottle Holder:

170mm

Nominal Weight:

210g (Bottle Holder & 500ml bottle)

Recommended Storage Conditions:

Dry and ambient temperature

Notes:

- 1 - Supplied with 1 x 500ml bottle (8126H-500).
- 2 - The bottle holder can be used with any round bottle from 25mm to 95mm outer diameter.

Spare Parts

The Bottle Holder will work with either of the following 500ml or 1000ml bottles

675 8126H500	Wide Mouth Bottle 500ml	Box of 48
675 8126H1000	Wide Mouth Bottle 1000ml	Box of 24

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 31-Jan-25