


1. Part Number: 675 5170A100

2. Photo:



### 3. Product Information:

Description	Utility Dipper 316L ss - 50ml
Material of Construction:	316L stainless steel
Nominal Sample Volume:	100ml
Surface Finish:	Better than 0.5 microns Ra
Outer Diameter of Sample Chamber:	32mm
Sample Chamber Dimensions:	



Overall Length:	490mm
Weight:	380g
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name:	Kuldeep Flora
Position:	Quality Manager
Date:	12-Dec-24