

**1. Part Number:** 675 A3661

**2. Photograph:**



### 3. Product Information:

|                                 |                                                                                       |
|---------------------------------|---------------------------------------------------------------------------------------|
| Description                     | Square PharmaScoop                                                                    |
| Nominal Volume:                 | 1500ml                                                                                |
| Material of Construction:       | 316L stainless steel                                                                  |
| Method of Construction:         | Crevice free, all welds ground and polished                                           |
| Surface Finish:                 | Better than 0.5 microns Ra                                                            |
| Overall Length:                 | 315 mm                                                                                |
| Body Length:                    | 190 mm                                                                                |
| Body Width:                     | 160 mm                                                                                |
| Body Height:                    | 48 mm                                                                                 |
| Nominal Weight:                 | 490 g                                                                                 |
| BSE/TSE statement:              | All polishing compounds used in the manufacture of this scoop are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature                                                           |

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Dec-24