

1. Part Number: 675 A376500

2. Photograph:



3. Product Information:

Description	Weighing Scoop 500ml
Nominal Volume:	500ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	320mm
Body Length:	205mm
Maximum Width of Scoop Body:	110mm
Minimum Width of Scoop Body:	65mm
Nominal Weight:	347g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Dec-24