

1. Part Number: 675 8015A170S

2. Photograph:



3. Product Information:

Description:	SteriWare Micro-Spatula - 170mm - Sterile
Use:	Single use, disposable powder spatula
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white • FDA acceptable CFR 177.1640 • EU Regulations 10/2011 compliant • EC Regulations 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Spatula Length:	175mm
Spatula Blade Length:	40mm
Blade Width:	13mm
Nominal Weight:	5g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45kGy)
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested