

1. Part Number: 675 8015H170S

2. Photograph:



3. Product Information:

Description:	SteriWare Micro Spatula HDPE - Sterile
Use:	Single use, disposable spatula
Material of Construction:	High Density Polyethylene (HDPE) - FDA acceptable CFR 177.1520 - EU Regulations 10/2011 compliant - EC Regulations 1935/2004 compliant - BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	2ml
Overall Length:	175mm
Spatula Blade Length:	40mm
Blade Width:	13mm
Nominal Weight:	4g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Scoops per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested