


1. Part Number: 675 8015N170

2. Photograph:



3. Product Information:

Description:	SteriWare Micro Spatula 0.5ml
Use:	Single use, disposable spatula
Material of Construction:	High Impact Polystyrene (HIPS) - FDA acceptable CFR 177.1640 - EU Regulations 10/2011 compliant - EC Regulations 1935/2004 compliant - BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	0.5ml
Overall Length:	170mm
Spatula Blade Length:	40mm
Blade Width:	6mm
Nominal Weight:	4g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Scoops per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested