

1. Part Number: 675 8015R170S

2. Photograph:



3. Product Information:

Description:	SteriWare Micro Spatula 2ml RED - Sterile
Use:	Single use, powder spatula
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none">• Conforms to FDA 21 CFR 177.1640• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	2ml
Overall Length	175mm
Blade Length:	40mm
Blade Width:	13mm
Nominal Weight:	5g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (Minimum 25 kGy)
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested