

1. Part Number: 675 8045A30

2. Photograph:



3. Product Information:

Description:

PharmaScoop 30ml

Use:

Single use, disposable scoop

Material of Construction:

High Impact Polystyrene (HIPS) – rigid, white

- FDA acceptable CFR 177.1640
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Overall Length:

153mm

Handle Length:

101mm

Maximum Width:

40mm

Nominal Weight:

10g

Moulding Environment:

Class 100,000 Cleanroom

Primary Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

N/A

Number of Scoops per Box:

100

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested