

1. Part Number: 675 8045B1000

2. Photograph:



3. Product Information:

Description:	SteriWare® PharmaScoop – 1000ml
Use:	Single use, disposable scoop
Material of Construction:	Polypropylene (PP) – virgin – blue • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	1000ml
Overall Length:	324mm
Handle Length:	149mm
Maximum Width:	150mm
Nominal Weight:	169g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	25
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested