

1. Part Number: 675 8045B125

2. Photograph:



3. Product Information:

Description:	PharmaScoop – 125ml Single use, disposable scoop Polypropylene
Use:	
Material of Construction:	(PP) – virgin – blue <ul style="list-style-type: none">• Conforms to FDA CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	125ml
Overall Length:	196mm
Handle Length:	104mm
Maximum Width:	61mm
Nominal Weight:	29g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested