

**1. Part Number:** 675 8045H250

**2. Photograph:**



### 3. Product Information:

Description:	PharmaScoop HDPE 250ml Single use,
Use:	disposable scoop
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"><li>• FDA acceptable CFR 177.1520</li><li>• EU Regulations 10/2011 compliant</li><li>• EC regulation 1935/2004 compliant</li><li>• BSE/TSE free</li></ul>
Method of Construction:	Single piece, injection moulded
Nominal Scoop Volume:	250ml
Overall Length:	209mm
Handle Length:	104mm
Maximum Width:	81mm
Nominal Weight:	38g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Scoops per Box:	85
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	The product has not been ATEX tested