

1. Part Number: 675 8045P1000

2. Photograph:



3. Product Information:

Description:

PharmaScoop PP 1000ml Single use,

Use:

disposable scoop Polypropylene –

Material of Construction:

rigid, opaque

- FDA acceptable CFR 177.1520
- EU Regulations 10/2011 compliant
- EC regulation 1935/2004 compliant
- BSE/TSE free

Method of Construction:

Single piece, injection moulded

Nominal Scoop Volume:

1000ml

Overall Length:

324mm

Handle Length:

149mm

Maximum Width:

150mm

Nominal Weight:

138g

Moulding Environment:

Class 100,000 Cleanroom

Primary Packaging Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

N/A

Number of Scoops per Box:

25

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

4. Change Control: