

1. Part Number: 675 8045P125

2. Photograph:



3. Product Information:

Description:	PharmaScoop PP 125ml Single use,
Use:	disposable scoop Polypropylene –
Material of Construction:	rigid, opaque • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Scoop Volume:	125ml
Overall Length:	196mm
Handle Length:	104mm
Maximum Width:	61mm
Nominal Weight:	27g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packaging Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Scoops per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested