


1. Part Number: 675 8045R2500

2. Photograph:



3. Product Information:

Description:	PharmaScoop RED 2500ml Single use,
Use:	disposable scoop
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, red • FDA acceptable CFR 177.1640 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	2500ml
Overall Length:	357mm
Handle Length:	147mm
Maximum Width:	166mm
Nominal Weight:	224g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packaging Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Scoops per Box:	15
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX rated