


1. Part Number: 675 8050H500

2. Photograph:



3. Product Information:

Description:

LiquiThief - 500 mm

Use:

Single use, disposable liquid sampler

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- Conforms to REACH regulation 1907/2006
- BSE/TSE free

O-ring Seal – virgin Silicone Rubber

- Conforms to FDA CFR 177.2600

Method of Construction:

Ultrasonic welding

Overall Length

Fully closed:

524mm

Fully extended:

885mm

Diameter of Outer Tube:

21mm

Diameter of hole in tip:

1.5mm

Width of Handle:

60mm

Nominal Weight of Sampler:

82g

Moulding Environment:

Class 100,000 Cleanroom

Assembly & Primary Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

N/A

Quantity per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

4. Change Control:

This is a single-use, disposable sampler. Repeated use will lead to failure.

5. Country of Origin:

6. Repeated Use:

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 2-Apr-25